




BIORESORBABLE FLAT FOOT IMPLANT




Surgical Technique
Pediatric
Treatment



NOVAGENIT®

INDICATIONS & CONTRAINDICATIONS

FLEXIBLE FLAT FOOT SURGICAL TREATMENT IN CHILDREN



Sandro Giannini M.D.

Villa Laura Clinic, Bologna, Italy

Patents:
Giannini US5360450
Giannini EP0560249B1

Indications Main indications:

- When the functional issue ceatly as-
sessed whith clinical examinaton, or in
uncertain cases, when it is confirmed by
gait analysis.
- Flat foot with tarsal fusion. In this case
after the fusion removal the implant
allows correction of the deformity
associated.
- Congential vertical talus. In this case
after soft tissue release, the device fills
the gap between talus and calca-
neous avoiding recuissence deformity.
- Adult flat foot. In this case the added
procedures (Achilles tendon leng-
thening and posterior tibial tendon
retention) are necessary.
- Posterior tibial tendon dysfunction
type 1° and 2°, associated with suture
repair of tendon transfer and Achilles
tendon lengthening when necessary.

CONTRAINDICATIONS

- Neurological flat foot
- Flat foot due to severe ligaments
laxity
- Adult flat foot with arthritis of the
hindfoot joints

ANAESTHESIA

Anaesthesia can be general,
peridural or peripheral, depending
upon the individual case. General
anaesthesia is usually preferred for
children.


Tourniquet:

The use of a tourniquet is optional.


OPERATIVE TECHNIQUE



The patient is placed in the supine position with the foot internally rotated. A 1cm incision over the sinus tarsi is performed.




The extensor retinaculum is opened to the cuboid bone using curved scissors. By turning the tip of the scissors upwards and pushing in a medial direction towards the internal malleolus, the tip of the scissors can be felt going into the deepest region of the sinus tarsi.




A 6mm rod is introduced in the same direction followed by 8 and 10mm rods until the correction is obtained.

The skin and the fibres of the retinaculum are opened with 2 small retractors to allow for the placement of the outer cylinder of the implant with a universal introducer.



The inner screw is inserted to open and stabilize the implant. The screw is tightened until the characteristic "squeaking" is heard from the material.



The retinaculum is sutured with N. 3-0 reabsorbable thread and the skin is sutured with another stitch.



Usually this is the only procedure necessary if the patient is still on the growth age.

ADDED SURGICAL PROCEDURES

After correcting the deformity by inserting the implant, the dorsiflexion of the foot is checked with the knee in extension. If the Achilles tendon is at 90°, the Achilles tendon is lengthened subcutaneously by two or three alternate hemisections, starting distally-laterally, then 3 cm above the previous one. The foot is forced in dorsal flexion to stretch the tendon until 10° dorsiflexion is achieved.



In case of navicular accessory or prominent painful navicular bone, or the interruption of the Meary's line at the naviculocuneiform with an angle greater than 10°, a medial procedure is performed with retention of the posterior tibialis.

With a 3cm incision over the navicular prominence, the navicular bone is revealed along with the posterior tibialis.

The periosteum is detached from the navicular bone and the posterior tibial tendon, maintaining the metatarsal expansion and

fibres directed towards the navicular bone.

After tangential resection of the navicular prominence, remove any navicular accessory.

The posterior tibialis is put under tension using a No. 2 reabsorbable stitch which passes through the dorsal periosteal flap.

Take the tendon of the posterior tibialis with a stitch, according Bunnell, and pass from plantar to dorsal through the spongy part of the navicular bone.

Added Surgical Procedures



By pulling the two ends of the thread, the tendon is advanced distally and fixed under the navicular bone. The suture is reinforced with another cross stitch.

The sheath of the posterior tibialis is sutured with No. 3-0 thread.

POSTOPERATIVE TREATMENT



When only the implant is inserted, it is preferable to use a walking boot for 2 weeks in order to reduce potential pain during walking and allow earlier return to normal activity. If other surgical procedures are associated, the recommended period of immobility with a boot is 5 weeks without weight bearing and 2 weeks with weight bearing. In both cases, when the boot is removed, normal footwear is worn and cycling and swimming are recommended.

CLINICAL CASES








Clinical aspect of an eleven-year-old boy before and after surgery and Rx before and after operation.

BEFORE



AFTER



MRI after 6 months shows the maintenance of the shape of the device, the device is breaking after 1.5 years, is partially reabsorbed after 3 years, and completely reabsorbed after 5 years.

REFERENCES

- **Benedetti, M., Ceccarelli, F., Berti, L., Luciani, D., Catani, F., Boschi, M., et al.** (2011). Diagnosis of flexible flatfoot in children: a systematic clinical approach. *Orthopedics*, p. 24(2):94.
- **Bertani, A.; Capello, A.; Catani, F.; Benedetti, M.G.; Simoncini, L.:** Flat foot Functional Evaluation Using Pattern Recognition Of Ground Reaction Data. *Clin. Biomech.*, 14: 484-493, 1999.
- **Bordelon, R.L.:** flexible Flatfoot. In: Mann RA, Coughlin MJ, eds. *Surgery of the Foot & Ankle*. Mosby S. Louis, 1993, pp. 717-756.
- **Carranza-Bencano, A., Duque-Gimeno, V., Gomez-Arroyo, J., & Zurita-Gutierrez, M.** (2000). Giannini's prosthesis in the treatment of the juvenile flatfoot. *Foot and Ankle Surg.*, p. 6:11-17.
- **Ceccarelli, F. B.** (2007). Le protesi di sottoastragalia - Subtalar point prosthesis. *G.I.O.T.*, p. 33:183-189.
- **De Rosa, G.P.:** Flexible Flatfoot. In: *Operative Foot Surgery*. Gould N, ed..WB Saunders, Philadelphia, 1994, pp. 834- 857.
- **Dockery, G.L.:** Symptomatic Juvenile Flatfoot Condition. *Surgical Treatment*. *J Foot Ankle Surg.*, 34: 135-145, 1995.
- **Giannini, S.; Ceccarelli, F.; Girolami, M.; Catani, F.; Benedetti, M.G.:** Kinematic and isokinetic evaluation of patients with flat foot. *Ital J Orthop Traumatol.*, 18: 249-260, 1992.
- **Giannini, S.; Girolami, M.; Ceccarelli, F.:** The surgical treatment of infantile flat foot. A new expanding endo-orthotic implant. *Ital J Orthop traumat.*, 11: 315-322, 1985
- **Giannini, S.:** Kenneth A. Johnson Memorial Lecture. Operative treatment of the flatfoot: why and how. *Foot Ankle Int.*, 19: 52-58, 1998
- **Giannini, S.; Ceccarelli F.:** The flexible flat foot. *Foot Ankle Clin.*, 3: 373-392, 1998
- **Giannini, S., Ceccarelli, F., Benedetti, M., Catani, F., & Faldini, C.** (2001). Surgical treatment of flexible flatfoot in children. *J Bone Joint Surgery*, p. 83A, 2, 2: 73-79.
- **Giannini, S., Ceccarelli, F., Vannini, F., & Baldi, E.** (2003). Operative treatment of flatfoot with talocalcaneal coalition. *Clin Orthop Relat Res*, p. 411:178-187.
- **Giannini, S., Faldini, C., Luciani, D., Cadossi, M., & Capra, P.** (2010). Piede piatto: endortesi riassorbibile – Flat foot: reabsorbable implant. *GIOT*, p. 36:278-284.
- **Giannini, S., Vannini, F., Faldini, C., Luciani, D., & Chiarello, E.** (2008). Endortesi nella patologia da anomala pronazione di sottoastragalia: vite ad espansione – Use of bioabsorbable implants for surgical treatment of flatfoot. *G.I.O.T.*, p. Supp.1:203-205.
- **Gutiérrez, P., & Herrera Lara, M.** (2005). Giannini Prosthesis for Flatfoot. *Foot & Ankle International*, p. 26:918.
- **Halabchi, F., Mazaheri, R., Mirshahi, M., & Abbasian, L.** (2013). Pediatric Flexible Flatfoot; Clinical Aspects and Algorithmic Approach. *Iran J Pediatr*, p. 23(3):247-260.
- **Harris, E., Vanore, J., & Thomas, J.** (2004). Diagnosis and treatment of pediatric flatfoot. *J Foot Ankle Surg*, p. 43(6):341-373.
- **Isikan, U.E.:** The values of talonavicular angles in patients with pes planus. *J. Foot Ankle Surg.*, 32: 514-516, 1993
- **Lee, M., Jo, H., & TS, B.** (2003). Analysis of initial fixation strength of press-fit fixation technique in anterior cruciate ligament reconstruction. A comparative study with titanium and bioabsorbable interference screw using porcine lower limb. *Knee Surg Sports Traumatol Arthrosc.*, p. 11(2): 91-98.
- **Root, M.L.; Orien, W.P.; Weed, J.H.:** Normal and abnormal function of the foot. *Clinical Biomechanics Corporation Publishers*, Los Angeles, 1977
- **Rose, G.K.; Welton, E.A.; Marshall, T.:** The Diagnosis Of Flat Foot In the Child. *J. Bone and Joint Surg.*, 67B: 71-78, 1985
- **Ruozzi, B., Belletti, D., Manfredini, G., Tonelli, M., Sena, P., Vandelli, M., et al.** (2013). Biodegradable device applied in flatfoot surgery: comparative study between clinical and technological aspects of removed screws. *Materials Science and Engineering*, p. 33: 1773-1782.
- **Vainionpää, S., Rokkanen, P., & Törmälä, P.** (1989). Surgical applications of biodegradable polymers in human tissues. *Prog. Polym. Sci.*, p. 14: 679.
- **Viladot, R.; Richera R.; Viladot A.:** Quince lecciones sobre patologia del pied. Barcelona, Ediciones Toray SA, 1989, pp. 69-93



IMPLANTS AND INSTRUMENTS



CODE IMPLANTS

DESCRIPTION

2100008

BFFI 8mm in sterile package

2100010

BFFI 10mm in sterile package

The flat foot implant is a Class III medical device

INSTRUMENTS (Not sterile Class I medical device)

FA00000

BFFI complete instrument set



FA00100

Cylindric spacer 6mm



FA00200

Cylindric spacer 8mm



FA00300

Cylindric spacer 10mm



FA00400

Implant introducer 8mm



FA00700

Implant introducer 10mm



FA00500

BFFI Screwdriver



Headquarters
viale Trento 115/117 38017
MEZZOLOMBARDO (TN) Italy
cod. fiscale e P. IVA 01949700221
Commercial Office

Via dell'Innovazione 9, 20032 Cormano, Milano, Italy
tel +39-02-6154 37404 - fax +39-02-6154 37417
www.novagenit.com
info@novagenit.com