



# FURLONG EVOLUTION

THE NATURAL SELECTION



# Contents

|   |    |
|---|----|
| The Furlong Legacy                          | 2  |
| Design Evolution                            | 4  |
| Furlong Evolution® Design Criteria          | 5  |
| Furlong Evolution® Proven Principles        | 6  |
| Furlong Evolution® Enhancements             | 7  |
| CSF <i>Plus</i> - High Performance Bearings | 8  |
| Pre-operative Planning                      | 10 |
| Femoral Neck Osteotomy                      | 11 |
| Acetabular Preparation                      | 12 |
| Acetabular Alignment                        | 13 |
| Cup Impaction                               | 14 |
| Screw Insertion                             | 15 |
| Femoral Preparation                         | 16 |
| Broaching the Femur                         | 17 |
| Trial Reduction                             | 18 |
| Liner Insertion                             | 19 |
| Stem Impaction                              | 20 |
| Final Reduction                             | 21 |
| Technical Data                              | 22 |
| Furlong Evolution Implant Range             | 23 |

“The preservation of favourable variations and the rejection of injurious variations, I call Natural Selection.”

**Charles Darwin**

On The Origin Of Species By Means of Natural Selection



**FURLONG EVOLUTION**


T H E   N A T U R A L  
S E L E C T I O N



# The Furlong® Legacy

By 1982 Mr Ronald Furlong F.R.C.S. had performed numerous revisions of cemented stems and had observed that the most common cause of implant failure was aseptic loosening. He concluded that to increase the longevity of a prosthesis, addressing this issue would be the number one priority. Collaborating with leading experts in biomechanics, engineering and biological implant fixation, he set out to develop an uncemented prosthesis with the principle objective of achieving:

- Absolute primary mechanical stability
- Rotational stability
- Long term physiological fixation



## **Rotational stability -**

rectangular cross section of the proximal body provides rotational stability under dynamic loading

## **Primary mechanical stability -**

In engineering the most efficient means of transferring load from one object to another is through the use of a cone. The cone prevents subsidence by creating 'hoop stresses' on the cortical bone

## **Secondary physiological**

**fixation -** Vacuum plasma sprayed Supravit® H-A.C. coating provides proven osseointegration and superior bond strength, resistance to delamination and resorption of the coating over time.

First implanted in September 1985, its long term clinical results have consistently performed and fulfilled Mr Furlongs criteria.

### Aseptic loosening – a thing of the past!



100%

16 – 19 years<sup>1</sup>

1. Survivorship of 38 cases in under 50 year olds. N.N. Shah *et al* **J Bone Joint Surg [Br]** 2009; 91-B:865-9



100%

15 – 21 years<sup>2</sup>


2. Survivorship of 331 consecutive cases. J.A.N Shepperd *et al* **J Bone Joint Surg [Br]** 2008; 90-B:27-30



100%

13 – 15 years<sup>3</sup>

3. Survivorship of 134 consecutive cases. A.A. Shetty *et al* **J Bone Joint Surg [Br]** 2005; 87-B:1050-4



99.73%

19 years<sup>4</sup>

4. Survivorship in 2,212 cases. J.M. Buchanan, Sunderland Royal Hospital  
Data presented at BOA 2007,  
Manchester, 26 - 28 September 2007



FURLONG EVOLUTION



## Design Evolution

---

To be successful it was acknowledged that a new design of uncemented femoral prosthesis must build on the strengths of the Furlong® H-A.C. with adaptations which optimise ease of use without compromising on the ability to achieve secure immediate and long term fixation and excellent clinical outcomes.

The Furlong® has evolved to even more closely match the needs of;

Patient

Surgeon

Theatre staff

Theatre management

Procurement

Sterile services





## Furlong Evolution® Design Criteria

---

- Maintain clinical results - preserve the principles set out by Mr Furlong to achieve continuing clinically successful outcomes
- Optimise surgical technique - to achieve accurate implantation with predictable and reproducible results
- Perfect the instrumentation - Simple, concise, easy to use



FURLONG EVOLUTION



# Furlong® Evolution Proven Principles

---



**12/14 macro taper** - unchanged on the Furlong® since 1989, the trunnion design allows the use of large head articulation without the risk of 'trunionosis'.

**Rectangular body cross section** - for rotational Stability.



**Cone transition between proximal body and cylindrical distal stem** - proven primary mechanical stability.

**Vacuum plasma sprayed Supravit® Hydroxyapatite ceramic coating** - proven long term fixation.





# Furlong Evolution® Enhancements

**Reduced lateral shoulder** - allows for easier insertion and optimum 'anatomical' fit.

**Shorter stem design** - preserves bone stock, facilitates contemporary surgical approaches and reduces soft tissue disruption.

**Extensive range of stem sizes** - 10 sizes from 8-17mm in 1mm increments, consistently increasing in both the antero-posterior and medio-lateral plane by 1mm per size, to provide optimum metaphyseal fit.

**Optimised neck geometry** - increasing the range of motion.



**Supravit® Zoned coated** - push-out tests have shown an increased shear strength at the H-A.C. to bone interface a function of surface roughness. Furlong Evolution® combines roughened coating proximally for improved 'scratch fit', while retaining the smoother Supravit® on the distal stem for easier insertion.

**Anatomic neck shaft angles** - all stem sizes available in two neck shaft angles 126° and 133° and two offsets, provides maximum versatility to restore joint biomechanics and soft tissue tension.

**Constant neck length and location on all stem sizes** - giving more predictable offsets

**Collared and non-collared** - all stem sizes are available with or without a collar

**Constant radius of the medial curve in all sizes** - conforms to the medial calcar curve enabling smooth insertion with predictable positioning'.



**FURLONG EVOLUTION**



## CSF *Plus* - High Performance Bearings

### BioloX Delta Ceramic



The Furlong® H-A.C. CSF *Plus* cup shell features Supravit® Zoned coating for secure initial fixation. The acetabular shell is designed with an 18° internal taper providing excellent fixation of the liner in the shell, as well as allowing easy removal should revision of the components ever be required. 36 mm internal diameter liners can be used from the 50mm outside diameter shell size upwards.

### CLP75® X-Link polyethylene



The liner options for the Furlong® H-A.C. CSF *Plus* polyethylene bearing cups are available in CLP75®, JRI's proprietary highly crosslinked polyethylene. To maximise polyethylene thickness each liner fits one shell size only. The 5° internal taper of the shell avoids the use of a potentially destructive locking mechanism yet allows for secure seating of the liner ensuring that the coupling is free of micro-movement, greatly reducing the risk of backside wear\*. 10° hooded liners are available for use with all shell sizes should additional stability be required. Where adjunctive screw fixation is not used, screw-hole plugs are provided to seal the cup against polyethylene debris migration.

\*Data presented at Efort 2003 (Helsinki) -The 6 to 10 year results of the CSF hydroxyapatite coated acetabular cup : Andrew F.M. McKee, M.D.George, R. Hussein, J.A.N. Shepperd : Conquest Hospital, Hastings, England

## TriboFit™ Polycarbonate Urethane (PCU)


The TriboFit™ Acetabular Buffer is manufactured from a hydrophilic material which helps to maintain a full fluid film layer between the buffer and the femoral head providing microelastohydrodynamic lubrication. In this situation, as in the natural joint there is no physical contact between the bearing surfaces and essentially no wear occurs.



- Low wear
- Cartilage-like material
- Resilient, soft, pliable bearing surface
- Shock absorbing
- Enhanced lubrication
- Decreased friction
- Increased stability
- Improved range of motion

‘The PCU liner showed excellent wear characteristics in terms of its low and steady volumetric wear rate (5.8-7.7 mm<sup>3</sup>/Mc) and low particle generation rate (2-3 × 10<sup>6</sup> particles/Mc). The latter is 5-6 orders of magnitude lower than that of highly cross-linked polyethylene and 6-8 orders of magnitude lower than that of metal-on-metal bearings.’

Long-term evaluation of a compliant cushion form acetabular bearing for hip joint replacement: a 20 million cycles wear simulation. Elsner JJ et al. Orthop Res. 2011 Dec;29(12):1859-66. doi: 10.1002/jor.21471. Epub 2011 May 31



Pre-operative planning, including X-ray templating, is strongly recommended for leg length and offset assessment and to provide guidance for accurate bone preparation and appropriate implant selection. A bi-lateral AP X-ray of the proximal femur and pelvis will aid in leg length and offset assessment and management.

Acetate radiographic templates for the Furlong Evolution® and CSF *Plus* are available in 115% magnification.

Digital templates are available in OrthoView™, mediCAD®, SECTRA® and TraumaCad® for use with PACS. If PACS digital templating is to be used, ensure the X-rays have been calibrated to the correct magnification using a reliable marker.


Caution: Pre-operative templating is intended for estimating purposes only. Final component size is determined intraoperatively.





# Femoral Neck Osteotomy

---



After dislocation, resect the femoral neck. The femoral neck resection guide can be used to aid alignment and inclination of the cut. The level of the femoral neck cut should be performed in conjunction with the pre-operative planning. The line of the femoral resection represents  $47^\circ$  from the long axis of the femur.

Stems are available in two neck shaft angles,  $126^\circ$  And  $133^\circ$ . Align the resection guide with the femoral axis and determine which neck shaft angle matches the patients anatomy.





## Acetabular Preparation

---



Image 2



Image 3

The acetabulum is exposed and reamed, starting with the smallest reamer (44mm) and increasing in 2mm increments until a bed of circumferential bleeding subchondral bone is achieved. The final reamer size determines the size of the definitive acetabular implant (Image 2).

The trial cup, corresponding to the final reamer, is screwed onto the introducer handle. The trial cup is introduced into the acetabulum to check for diameter and depth. The height of the trial cup is the same as the definitive implant, whereas the diameter of the trial mimics the reamer and NOT the definitive implant (Image 3).

*Note: The CSF Plus cups are hemispherical and are 2.5% larger in diameter at the equator than their corresponding reamers and test cups, plus the thickness of the H-A.C. coating (350µm). This facilitates a secure press fit when reamed line to line.*

In hard sclerotic bone it may be necessary to over ream the acetabular rim 1 or 2mm to aid insertion. Do not fully over-ream the socket.



# Acetabular Alignment



Image 4




Image 5




Image 6

The definitive acetabular implant is firmly screwed onto the impactor shaft. Align the cluster of 3 holes in the acetabular cup to the superior aspect of the prepared acetabulum (Image 4). The isolated holes should now line up with the ischium and the pubis.

The cup positioning guide is clipped onto the impactor shaft and secured.

The horizontal guide rod is screwed into the body of the guide in the corresponding left or right hole depending on which side is being operated on. The offset holes represent 20° of anteversion when aligned with the long axis of the body (Image 5). The vertical and horizontal rods, when perpendicular and parallel to the floor respectively, indicate an abduction angle of 45° (Image 6).





## Cup Impaction

---



Image 7



Image 8

The cup is presented to the prepared acetabulum ensuring the three clustered screw holes align with the ilium (Image 7). Impact the cup firmly ensuring the orientation of the cup is maintained. An audible change of tone and tactile feedback will indicate when the cup is correctly seated. The impactor handle can be removed and the seating of the cup checked through the apical hole and/or through the screw holes (Image 8).

Should the cup be difficult to seat or the bone be sclerotic, it may be necessary to open the aperture of the acetabular rim. Lightly reaming the acetabular rim with the next size reamer up will help to introduce and seat the implant.





## Screw Insertion



Image 9



Image 10



Image 11



Image 12

Additional adjunctive fixation can be achieved with the use of cancellous bone screws to supplement primary mechanical fixation.

A drill guide is placed carefully in the chosen screw hole, ensuring that it is fully seated. A pilot hole is made using the disposable drill bit which can be used with power or manually (Image 9).

The length of the screw is determined using the depth gauge (Image 10). Once selected the self tapping screws, which are available in 5mm increments from 15mm to 50mm, are held with the screw holding forceps and screwed into position with a 3.5mm hex screwdriver (Image 11). It is imperative that the head of the screw does not protrude from the screw hole, failure to do so could result in the screw head preventing the ceramic liner from seating fully (Image 12).

Screws can be angled through 30° to achieve optimum fixation with the cancellous bone.

The trial insert with the desired internal diameter is fitted to prevent any damage to the internal taper.





## Femoral Preparation

---



Image 13



Image 14

The box chisel attaches to the broach handle and is placed just medial to the resected cortical bone, within the boundary of the cancellous region of the medial calcar. An initial wedge of cancellous bone is resected using the box chisel whilst keeping it in line with the femur. If excessive version is required later in the procedure it is advisable to consider this when making the cut with the box chisel (Image 13).

The small tapered t-handled intramedullary reamer is used to open up the medullary canal. A larger reamer is available should the femoral canal be wide. Both reamers have a laser mark to indicate the minimum depth. The reamer is introduced until the line is at minimum, level with where the proximal shoulder of the prosthesis is intended to be as indicated in the pre-operative planning (Image 14).



## Broaching the Femur



Image 15



Image 16




Image 17

Attach the broach to the broach handle and prepare the femoral canal. Beginning at size 8mm, progress through the available sizes until the desired implant size is reached. To reduce the risk of damage to the calcar, care should be taken to gently and repeatedly hammer the broach in and out (Image 15). The proximal etched line indicates the 'stop' point of the broach. This line should be visible above the resection line medially (Image 17).





## Trial Reduction



Image 18



Image 19

Upon reaching the desired implant size detach the broach handle from the broach in order to facilitate use of plastic, colour coded trial necks.

Trial necks are available in 4 clearly labelled options (Image 19):

- 133deg standard offset (yellow)
- 133deg high offset (yellow)
- 126deg standard offset (orange)
- 126deg high offset (orange)

Attach a trial neck to the broach and the appropriate trial head as templated during pre-operative planning. A choice of four neck lengths are available: Short (-4mm), Medium (0mm), Long (+4mm) and X-Long (+8mm). It is vital to select head diameter size which matches the chosen acetabular cup/liner I.D. Size (Image 18). Reduce the hip and check for stability and range of motion.

*Note: The black laser line on the trial neck indicates that correct orientation when visible from above.*



## Liner Insertion




Image 20



Image 21



Image 22



Image 23



Image 24

Remove the trial head, neck, broach and trial liner, removal forceps are available to assist in removing the trial liner (Image 20).


The inner surface of the acetabular cup must be clean, dry and free from debris. The definitive acetabular liner is inserted into the cup by gently sliding the liner down the internal taper of the cup (Images 21 and 22). Ensure that the liner is uniformly flush against the whole rim by gliding your finger around the edge of the cup confirming the edge of the liner is even with the rim of the shell around the whole circumference (Image 23).

The liner is firmly tapped home using the trial cup introducer handle with the size of plastic cap that matches the internal diameter of the acetabular liner. After impaction re-check to ensure correct seating of the liner.





## Stem Impaction



The definitive femoral stem is implanted using the stem introducer. Two introducer options are available: threaded and hexagon ended. Both styles of stem introducer locate into the threaded/hexagon hole lateral to the shoulder of the definitive implant. It is advisable to insert the stem as far as possible by hand (within 10mm of the final position). Final impaction can be achieved with the aid of light tapping using a mallet. A change of tone can be used as a guide to determine the final position of the stem (Image 25).

- If using a collared stem option, the collar is designed to sit 1mm off the calcar. A change in tone will indicate when the stem is seated. Further impaction to attempt to seat the collar on the calcar is not necessary (Image 26).
- If using a non-collared stem final seating is achieved when the H-A.C. coating is in line with the neck resection line. Again the note may change to a harder tone before the desired position is achieved and no further impaction is necessary (Image 27).



## Final Reduction



Image 28



Image 29



Image 30



Image 31

Once the stem is fully impacted, the spigot protector is removed. A final trial reduction is performed to confirm joint stability and range of motion (Image 28).

The definitive femoral head with the desired neck length is placed carefully onto the clean spigot of the femoral stem using a twisting motion. Using a light tap, impact the definitive femoral head using the head impactor (Image 29), ensuring the surface of the head is not scratched or damaged in any way.

The hip is reduced (Image 31).





# Femoral Stem Removal

---



Image 32



Image 33



Image 34


Should the stem require removal, slide the slaphammer weight onto the thinner section of the introducer handle closest to the thread. This must be done before the handle is screwed onto the stem (image 32).

The introducer handle is screwed firmly into the stem (image 33), the stem should be removed with caution in a curved manner as not to damage the greater trochanter (image 34).





## Technical Data



**A** = 100mm (All stem sizes)

**B** = Stem size (Distal stem diameter)

**G** = Proximal body lateromedial width (Stem size (B) + 2mm)

|         | Offset (D) | Neck Length (E) | Neck Height (F) |
|---------|------------|-----------------|-----------------|
| 133°    | 36.9       | 32              | 31              |
| 133° HO | 41.9       | 35.7            | 31              |
| 126°    | 40.9       | 31.6            | 26.4            |
| 126° HO | 45.9       | 35.3            | 26.4            |

| Offset value with varying neck lengths |        |       |        |        |             |             |
|--|--------|-------|--------|--------|-------------|-------------|
|  | S (-4) | M (0) | L (+4) | L (+8) | XL +4 (+12) | XL +8 (+16) |
| 133°                                   | 34     | 36.9  | 39.9   | 42.8   | 45.7        | 48.6        |
| 133° HO (+5mm)                         | 39     | 41.9  | 44.9   | 47.8   | 50.7        | 53.6        |
| 126°                                   | 37.7   | 40.9  | 44.2   | 47.4   | 50.6        | 53.9        |
| 126° HO (+5mm)                         | 42.7   | 45.9  | 49.2   | 52.4   | 55.6        | 58.9        |

All measurements in mm unless otherwise stated



FURLONG EVOLUTION



# Furlong Evolution Implant Range



## Collared, 126° CCD Angle, Standard Offset

|           |         |
|-----------|---------|
| 426.01.08 | Size 8  |
| 426.01.09 | Size 9  |
| 426.01.10 | Size 10 |
| 426.01.11 | Size 11 |
| 426.01.12 | Size 12 |
| 426.01.13 | Size 13 |
| 426.01.14 | Size 14 |
| 426.01.15 | Size 15 |
| 426.01.16 | Size 16 |
| 426.01.17 | Size 17 |



## Collared, 126° CCD Angle, High Offset

|           |         |
|-----------|---------|
| 426.51.08 | Size 8  |
| 426.51.09 | Size 9  |
| 426.51.10 | Size 10 |
| 426.51.11 | Size 11 |
| 426.51.12 | Size 12 |
| 426.51.13 | Size 13 |
| 426.51.14 | Size 14 |
| 426.51.15 | Size 15 |
| 426.51.16 | Size 16 |
| 426.51.17 | Size 17 |



## Collared, 133° CCD Angle, Standard Offset

|           |         |
|-----------|---------|
| 433.01.08 | Size 8  |
| 433.01.09 | Size 9  |
| 433.01.10 | Size 10 |
| 433.01.11 | Size 11 |
| 433.01.12 | Size 12 |
| 433.01.13 | Size 13 |
| 433.01.14 | Size 14 |
| 433.01.15 | Size 15 |
| 433.01.16 | Size 16 |
| 433.01.17 | Size 17 |



## Collared, 133° CCD Angle, High Offset

|           |         |
|-----------|---------|
| 433.51.08 | Size 8  |
| 433.51.09 | Size 9  |
| 433.51.10 | Size 10 |
| 433.51.11 | Size 11 |
| 433.51.12 | Size 12 |
| 433.51.13 | Size 13 |
| 433.51.14 | Size 14 |
| 433.51.15 | Size 15 |
| 433.51.16 | Size 16 |
| 433.51.17 | Size 17 |



**Collarless, 126° CCD  
Angle, Standard Offset**

|           |         |
|-----------|---------|
| 426.02.08 | Size 8  |
| 426.02.09 | Size 9  |
| 426.02.10 | Size 10 |
| 426.02.11 | Size 11 |
| 426.02.12 | Size 12 |
| 426.02.13 | Size 13 |
| 426.02.14 | Size 14 |
| 426.02.15 | Size 15 |
| 426.02.16 | Size 16 |
| 426.02.17 | Size 17 |



**Collarless, 126° CCD  
Angle, High Offset**

|           |         |
|-----------|---------|
| 426.52.08 | Size 8  |
| 426.52.09 | Size 9  |
| 426.52.10 | Size 10 |
| 426.52.11 | Size 11 |
| 426.52.12 | Size 12 |
| 426.52.13 | Size 13 |
| 426.52.14 | Size 14 |
| 426.52.15 | Size 15 |
| 426.52.16 | Size 16 |
| 426.52.17 | Size 17 |



**Collarless, 133° CCD  
Angle, Standard Offset**

|           |         |
|-----------|---------|
| 433.02.08 | Size 8  |
| 433.02.09 | Size 9  |
| 433.02.10 | Size 10 |
| 433.02.11 | Size 11 |
| 433.02.12 | Size 12 |
| 433.02.13 | Size 13 |
| 433.02.14 | Size 14 |
| 433.02.15 | Size 15 |
| 433.02.16 | Size 16 |
| 433.02.17 | Size 17 |



**Collarless, 133° CCD  
Angle, High Offset**

|           |         |
|-----------|---------|
| 433.52.08 | Size 8  |
| 433.52.09 | Size 9  |
| 433.52.10 | Size 10 |
| 433.52.11 | Size 11 |
| 433.52.12 | Size 12 |
| 433.52.13 | Size 13 |
| 433.52.14 | Size 14 |
| 433.52.15 | Size 15 |
| 433.52.16 | Size 16 |
| 433.52.17 | Size 17 |





# Polyethylene CLP75® Bearing Cup Range

H-A.C. CSF *Plus* - Liner/head  
size combination chart  
CLP75® Polyethylene



Standard Liners (I/D)

| Cup/Reamer<br>size (O/D) | 28mm      | 32mm      | 36mm      |
|--------------------------|-----------|-----------|-----------|
| 44mm<br>175.44.07        | 150.44.28 |           |           |
| 46mm<br>175.46.07        | 150.46.28 |           |           |
| 48mm<br>175.48.07        | 150.48.28 | 150.48.32 |           |
| 50mm<br>175.50.07        | 150.50.28 | 150.50.32 |           |
| 52mm<br>175.52.07        | 150.52.28 | 150.52.32 |           |
| 54mm<br>175.54.07        | 150.54.28 | 150.54.32 |           |
| 56mm<br>175.56.07        | 150.56.28 | 150.56.32 | 150.56.36 |
| 58mm<br>175.58.07        | 150.58.28 | 150.58.32 | 150.58.36 |
| 60mm<br>175.60.07        | 150.60.28 | 150.60.32 | 150.60.36 |
| 62mm<br>175.62.07        | 150.62.28 | 150.62.32 | 150.62.36 |



|        |          |          |          |
|--------|----------|----------|----------|
| Short  | 47.28.10 | 47.32.10 | 47.36.10 |
| Medium | 47.28.20 | 47.32.20 | 47.36.20 |
| Long   | 47.28.30 | 47.32.30 | 47.36.30 |
| X Long | 47.28.40 | 47.32.40 | 47.36.40 |

#### Hooded Liners (I/D)

| 28mm      | 32mm      | 36mm      |
|-----------|-----------|-----------|
| 151.44.28 |           |           |
| 151.46.28 |           |           |
| 151.48.28 | 151.48.32 |           |
| 151.50.28 | 151.50.32 |           |
| 151.52.28 | 151.52.32 |           |
| 151.54.28 | 151.54.32 |           |
| 151.56.28 | 151.56.32 | 151.56.36 |
| 151.58.28 | 151.58.32 | 151.58.36 |
| 151.60.28 | 151.60.32 | 151.60.36 |
| 151.62.28 | 151.62.32 | 151.62.36 |

|        |           |           |           |
|--------|-----------|-----------|-----------|
| Short  | 192.28.10 | 192.32.10 | 192.36.10 |
| Medium | 192.28.20 | 192.32.20 | 192.36.20 |
| Long   | 192.28.30 | 192.32.30 | 192.36.30 |
| X Long |           |           | 192.36.40 |





## CSF *Plus* Ceramic Bearing Cup Range

| Cup/Reamer size (O/D) | 28mm      | 32mm      | 36mm      | 40mm      |
|-----------------------|-----------|-----------|-----------|-----------|
| 44mm<br>170.44.07     | 150.44.28 |           |           |           |
| 46mm<br>170.46.07     |           | 140.46.32 |           |           |
| 48mm<br>170.48.07     |           |           |           |           |
| 50mm<br>170.50.07     |           | 140.50.32 | 140.50.36 |           |
| 52mm<br>170.52.07     |           |           |           |           |
| 54mm<br>170.54.07     |           | 140.54.32 | 140.54.36 | 140.54.40 |
| 56mm<br>170.56.07     |           |           |           |           |
| 58mm<br>170.58.07     |           |           |           |           |
| 60mm<br>170.60.07     |           |           |           |           |
| 62mm<br>170.62.07     |           |           |           |           |



H-A.C. CSF *Plus* - Liner/head size combination chart



Bilox® *delta* Ceramic





## TriboFit™ Implant Range



### Acetabular shell, Supravit® Zoned coated Cobalt Chrome

|           |          |
|-----------|----------|
| 1007 3046 | 46mm O/D |
| 1007 3048 | 48mm O/D |
| 1007 3050 | 50mm O/D |
| 1007 3052 | 52mm O/D |
| 1007 3054 | 54mm O/D |
| 1007 3056 | 56mm O/D |
| 1007 3058 | 58mm O/D |
| 1007 3060 | 60mm O/D |
| 1007 3062 | 62mm O/D |



### Acetabular Buffer Polycarbonate Urethane

|           |                    |
|-----------|--------------------|
| 1001 0040 | 40mm O/D, 34mm I/D |
| 1001 0042 | 42mm O/D, 36mm I/D |
| 1001 0044 | 44mm O/D, 38mm I/D |
| 1001 0046 | 46mm O/D, 40mm I/D |
| 1001 0048 | 48mm O/D, 42mm I/D |
| 1001 0050 | 50mm O/D, 44mm I/D |
| 1001 0052 | 52mm O/D, 46mm I/D |
| 1001 0054 | 54mm O/D, 48mm I/D |
| 1001 0056 | 56mm O/D, 50mm I/D |

### Femoral Heads Cobalt Chrome

|           |                      |
|-----------|----------------------|
| 1002 034S | 34mm O/D Short Neck  |
| 1002 034M | 34mm O/D Medium Neck |
| 1002 034L | 34mm O/D Long Neck   |
| 1002 036S | 36mm O/D Short Neck  |
| 1002 036M | 36mm O/D Medium Neck |
| 1002 036L | 36mm O/D Long Neck   |





# TriboFit™ Implant Range



## Modular Femoral Heads Cobalt Chrome

|           |          |
|-----------|----------|
| 1002 0038 | 38mm O/D |
| 1002 0040 | 40mm O/D |
| 1002 0042 | 42mm O/D |
| 1002 0044 | 44mm O/D |
| 1002 0046 | 46mm O/D |
| 1002 0048 | 48mm O/D |
| 1002 0050 | 50mm O/D |



## Femoral Neck Adaptors Titanium

For 38-44mm O/D Femoral Heads

|           |             |
|-----------|-------------|
| 1003 001S | Short Neck  |
| 1003 001M | Medium Neck |
| 1003 001L | Long Neck   |

## Femoral Neck Adaptors Titanium

For 46-50mm O/D Femoral Heads

|           |             |
|-----------|-------------|
| 1003 002S | Short Neck  |
| 1003 002M | Medium Neck |
| 1003 002L | Long Neck   |





# Furlong Evolution Instrumentation

---

| Part Number | Description                            |
|-------------|--|
| 10-99-47    | Head Impactor handle                   |
| 10-99-38    | Head Impactor                          |
| 10-00-46    | Cushion, impactor                      |
| 70-28-10    | Trial Heads - 28S                      |
| 70-28-20    | Trial Heads - 28M                      |
| 70-28-30    | Trial Heads - 28L                      |
| 70-28-40    | Trial Heads - 28XL                     |
| 70-32-10    | Trial Heads - 32S                      |
| 70-32-20    | Trial Heads - 32M                      |
| 70-32-30    | Trial Heads - 32L                      |
| 70-32-40    | Trial Heads - 32XL                     |
| 70-36-10    | Trial Heads - 36S                      |
| 70-36-20    | Trial Heads - 36M                      |
| 70-36-30    | Trial Heads - 36L                      |
| 70-36-40    | Trial Heads - 36XL                     |
| 70-40-10    | Trial Heads - 40S                      |
| 70-40-20    | Trial Heads - 40M                      |
| 70-40-30    | Trial Heads - 40L                      |
| 70-40-40    | Trial Heads - 40XL                     |
| 86-99-14    | Curved stem Introducer                 |
| 450-99-01   | Instrument tray                        |
| 450-04-00   | Offset Broach Handle OPTION            |
| 450-04-01   | Straight Broach Handle OPTION          |
| 450-04-08   | Tommy Bar                              |
| 450-14-00   | 'Namba' hip slide                      |
| 450-06-08   | 8mm broach                             |
| 450-06-09   | 9mm broach                             |
| 450-06-10   | 10mm broach                            |
| 450-06-11   | 11mm broach                            |
| 450-06-12   | 12mm broach                            |
| 450-06-13   | 13mm broach                            |
| 450-06-14   | 14mm broach                            |
| 450-06-15   | 15mm broach                            |
| 450-06-16   | 16mm broach                            |
| 450-05-17   | 17mm broach                            |
| 450-99-1    | Instrument tray                        |
| 450-03-00   | Tapered Reamer Small                   |
| 450-03-01   | Tapered Reamer Large                   |
| 450-02-01   | Box Chisel (7.5mm)                     |
| 450-07-00   | Stem Inserter/Extractor                |
| 450-01-00   | Femoral Resection guide                |
| 450-10-00   | Slap Hammer                            |
| 452-00-26   | Trial Neck - 126° (orange)             |
| 452-00-33   | Trial Neck - 133° (yellow)             |
| 452-05-26   | Trial Neck - High Offset 126° (orange) |
| 452-05-33   | Trial Neck - High Offset 133° (yellow) |





# CSF *Plus* Instrumentation

---

## Instrumentation for Furlong® H-A.C. CSF *Plus* Acetabular Cup

| Part Number | Description                            |
|-------------|--|
| 10.00.26    | Hudson adaptor                         |
| 10.02.00    | Acetabular Reamer Shaft                |
| 10.04.44    | Acetabular Reamers 44mm O/D            |
| 10.04.46    | Acetabular Reamers 46mm O/D            |
| 10.04.47    | Acetabular Reamers 47mm O/D            |
| 10.04.48    | Acetabular Reamers 48mm O/D            |
| 10.04.50    | Acetabular Reamers 50mm O/D            |
| 10.04.52    | Acetabular Reamers 52mm O/D            |
| 10.04.54    | Acetabular Reamers 54mm O/D            |
| 10.04.56    | Acetabular Reamers 56mm O/D            |
| 10.04.58    | Acetabular Reamers 58mm O/D            |
| 10.04.60    | Acetabular Reamers 60mm O/D            |
| 10.04.62    | Acetabular Reamers 62mm O/D            |
| 11.10.44    | Test Cups 44mm                         |
| 11.10.46    | Test Cups 46mm                         |
| 11.10.48    | Test Cups 48mm                         |
| 11.10.50    | Test Cups 50mm                         |
| 11.10.52    | Test Cups 52mm                         |
| 11.10.54    | Test Cups 54mm                         |
| 11.10.56    | Test Cups 56mm                         |
| 11.10.58    | Test Cups 58mm                         |
| 11.10.60    | Test Cups 60mm                         |
| 11.10.62    | Test Cups 62mm                         |
| 10.99.26    | Cup introducer                         |
| 10.28.53    | Plastic Caps for Liner Impactor 28mm   |
| 10.32.53    | Plastic Caps for Liner Impactor 32mm   |
| 10.36.53    | Plastic Caps for Liner Impactor 36mm   |
| 10.40.53    | Plastic Caps for Liner Impactor 40mm   |
| 50.00.30    | Screw holding forceps                  |
| 50.00.37    | Flat forceps for removing trial liners |
| 50.00.60    | Flexible aw                            |
| 50.91.35    | 3.5mm flexible screwdriver             |
| 50.92.35    | 3.5mm hexagon screwdriver              |
| 50.91.00    | Ratchet handle                         |
| 55.00.00    | Screw depth gauge                      |
| 81.00.44    | Instrument tray                        |
| 86.11.03    | Cup introducer                         |
| 86.00.02    | Flexible drill shaft                   |
| 86.00.07    | Flexible drill guide                   |
| 86.11.09    | Cup positioning guide                  |



# Trial Polyethylene Liners for CSF *Plus*

## Trial Liners for CSF *Plus* Ceramic Liners Only

| Part Number | Description   |
|-------------|---|
| 240.44.28   | Trial Liners 28mm I/D, fits 44mm O/D                  |
| 240.46.32   | Trial Liners 32mm I/D, fits 46 & 48mm O/D             |
| 240.50.36   | Trial Liners 36mm I/D, fits 50 & 52mm O/D             |
| 240.54.40   | Trial Liners 40mm I/D, fits 54, 56, 58, 60 & 62mm O/D |
| 81.00.45    | Instrument Tray                                       |

## Instrumentation for Furlong® H-A.C. CSF *Plus* UHMWPE Bearing Cup

| Part Number | Description                                   |
|-------------|---|
| 250.46.28   | Standard Trial Liners 28mm I/D, fits 46mm O/D |
| 250.48.28   | Standard Trial Liners 28mm I/D, fits 48mm O/D |
| 250.50.28   | Standard Trial Liners 28mm I/D, fits 50mm O/D |
| 250.52.28   | Standard Trial Liners 28mm I/D, fits 52mm O/D |
| 250.54.28   | Standard Trial Liners 28mm I/D, fits 54mm O/D |
| 250.56.28   | Standard Trial Liners 28mm I/D, fits 56mm O/D |
| 250.58.28   | Standard Trial Liners 28mm I/D, fits 58mm O/D |
| 250.60.28   | Standard Trial Liners 28mm I/D, fits 60mm O/D |
| 250.62.28   | Standard Trial Liners 28mm I/D, fits 62mm O/D |
| 250.48.32   | Standard Trial Liners 32mm I/D, fits 48mm O/D |
| 250.50.32   | Standard Trial Liners 32mm I/D, fits 50mm O/D |
| 250.52.32   | Standard Trial Liners 32mm I/D, fits 52mm O/D |
| 250.54.32   | Standard Trial Liners 32mm I/D, fits 54mm O/D |
| 250.56.32   | Standard Trial Liners 32mm I/D, fits 56mm O/D |
| 250.58.32   | Standard Trial Liners 32mm I/D, fits 58mm O/D |
| 250.60.32   | Standard Trial Liners 32mm I/D, fits 60mm O/D |
| 250.62.32   | Standard Trial Liners 32mm I/D, fits 62mm O/D |

## Instrumentation for Furlong® H-A.C. CSF *Plus* UHMWPE Bearing Cup

| Part Number | Description                                   |
|-------------|---|
| 250.48.32   | Standard Trial Liners 32mm I/D, fits 48mm O/D |
| 250.50.32   | Standard Trial Liners 32mm I/D, fits 50mm O/D |
| 250.52.32   | Standard Trial Liners 32mm I/D, fits 52mm O/D |
| 250.54.32   | Standard Trial Liners 32mm I/D, fits 54mm O/D |
| 250.56.32   | Standard Trial Liners 32mm I/D, fits 56mm O/D |
| 250.58.32   | Standard Trial Liners 32mm I/D, fits 58mm O/D |
| 250.60.32   | Standard Trial Liners 32mm I/D, fits 60mm O/D |
| 250.62.32   | Standard Trial Liners 32mm I/D, fits 62mm O/D |
| 250.56.36   | Standard Trial Liners 36mm I/D, fits 56mm O/D |
| 250.58.36   | Standard Trial Liners 36mm I/D, fits 58mm O/D |
| 250.60.36   | Standard Trial Liners 36mm I/D, fits 60mm O/D |
| 250.62.36   | Standard Trial Liners 36mm I/D, fits 62mm O/D |
| 250.60.40   | Standard Trial Liners 40mm I/D, fits 60mm O/D |
| 250.62.40   | Standard Trial Liners 40mm I/D, fits 62mm O/D |





# TriboFit™ Instrumentation

---

| Part Number | Description                                |
|-------------|--|
| 1101 0001   | TriboFit General Instruments Base          |
| 1101 9999   | TriboFit General Instruments Case Lid      |
| 1151 034S   | 34mm Modular Trial Head - Short            |
| 1151 034M   | 34mm Modular Trial Head - Med              |
| 1151 034L   | 34mm Modular Trial Head - Long             |
| 1151 036S   | 36mm Modular Trial Head - Short            |
| 1151 036M   | 36mm Modular Trial Head - Med              |
| 1151 036L   | 36mm Modular Trial Head - Long             |
| 1151 0038   | 38mm Modular Trial Head                    |
| 1151 0040   | 40mm Modular Trial Head                    |
| 1151 0042   | 42mm Modular Trial Head                    |
| 1151 0044   | 44mm Modular Trial Head                    |
| 1151 0046   | 46mm Modular Trial Head                    |
| 1151 0048   | 48mm Modular Trial Head                    |
| 1151 0050   | 50mm Modular Trial Head                    |
| 1151 0099   | Femoral Head Impactor                      |
| 1151 0199   | Taper Insert Removal Tool                  |
| 1154 0000   | Spin Test Twister Handle                   |
| 1152 0015   | Threaded Neck Adapter Trial - Short (-4mm) |
| 1152 0016   | Threaded Neck Adapter Trial - Med (0mm)    |
| 1152 0017   | Threaded Neck Adapter Trial - Long (+4mm)  |
| 1154 0001   | Twister for Small Adaptor Heads (38-44mm)  |
| 1154 0002   | Twister for Medium Adaptor Heads (46-50mm) |
| 1154 0003   | Twister for Large Adaptor Heads (52-58mm)  |
| 1154 0003   | Twister for Large Adaptor Heads (52-58mm)  |
| 1101 0002   | TriboFit® General Instruments Tray         |
| 1153 0040   | 40mm Buffer Placer                         |
| 1153 0042   | 42mm Buffer Placer                         |
| 1153 0044   | 44mm Buffer Placer                         |
| 1153 0046   | 46mm Buffer Placer                         |
| 1153 0048   | 48mm Buffer Placer                         |
| 1153 0050   | 50mm Buffer Placer                         |
| 1153 0052   | 52mm Buffer Placer                         |
| 1153 0054   | 54mm Buffer Placer                         |
| 1153 0056   | 56mm Buffer Placer                         |
| 1101 0000** | TriboFit® General Instruments Case         |



## TriboFit® Trial Head Set

| Part Number | Description                   |
|-------------|-------------------------------|
| 91.16.34    | Tribofit Trial Head 34mm      |
| 91.16.36    | Tribofit Trial Head 36mm      |
| 91.16.38    | Tribofit Trial Head 38mm      |
| 91.16.40    | Tribofit Trial Head 40mm      |
| 91.16.42    | Tribofit Trial Head 42mm      |
| 91.16.44    | Tribofit Trial Head 44mm      |
| 91.16.46    | Tribofit Trial Head 46mm      |
| 91.16.48    | Tribofit Trial Head 48mm      |
| 91.16.50    | Tribofit Trial Head 50mm      |
| 91.15.10    | Tribofit Head Neck Adaptor -4 |
| 91.15.20    | Tribofit Head Neck Adaptor 0  |
| 91.15.30    | Tribofit Head Neck Adaptor +4 |

## TriboFit® Shell Impactor Set

| Part Number | Description                           |
|-------------|---------------------------------------|
| 1156-4046   | 46mm Acetabular Shell Impactor Collar |
| 1156-4048   | 48mm Acetabular Shell Impactor Collar |
| 1156-4050   | 50mm Acetabular Shell Impactor Collar |
| 1156-4052   | 52mm Acetabular Shell Impactor Collar |
| 1156-4054   | 54mm Acetabular Shell Impactor Collar |
| 1156-4056   | 56mm Acetabular Shell Impactor Collar |
| 1156-4058   | 58mm Acetabular Shell Impactor Collar |
| 1156-4060   | 60mm Acetabular Shell Impactor Collar |
| 1156-4062   | 62mm Acetabular Shell Impactor Collar |
| 1156-1099   | Acetabular Shell Impactor Handle      |
| 50.04.08    | Multi Size Allen Key                  |





# Acknowledgements

---

JRI Orthopaedics would like to thank the clinicians who have contributed to the design of the Furlong Evolution®, whose input and assessment have driven the development of this prosthesis.

## Professor Justin Cobb

Consultant Orthopaedic Surgeon  
Imperial College London

## Mr Johan Witt

Consultant Orthopaedic Surgeon  
University College London

## Mr Matthew Burwell

Consultant Orthopaedic Surgeon  
Royal United Hospital Bath

## Mr Jonathan Walczack

Consultant Orthopaedic Surgeon  
South London Healthcare Trust

## Mr Raghu Raman

Consultant Orthopaedic Surgeon  
Hull Royal Infirmary

## Mr Kerry Acton

Consultant Orthopaedic Surgeon  
Mount Alvernia Hospital

## Dr Gunnar Flivik

Associate Professor  
Skåne University Hospital, Sweden

## Mr Adeel Aqil

Research Fellow to Prof Cobb  
Imperial College London

Disclaimer: No royalties have been paid or will be paid for the development of this prosthesis, JRI Orthopaedics is wholly owned by Orthopaedic Research UK (ORUK) to whom it gift aids its profit, after operational costs. ORUK independently funds orthopaedic research for the advancement in knowledge for the prevention and treatment of bone and joint disease.

---

### Mr Tony Ward

Consultant Orthopaedic Surgeon  
Frenchay Hospital

### Mr Joyti Saksena

Consultant Orthopaedic Surgeon  
The Whittington Hospital

### Mr Lee Taylor

Consultant Orthopaedic Surgeon  
St Richard's Hospital

### Mr Ian Bacarese-Hamilton

Consultant Orthopaedic Surgeon  
The Whittington Hospital

### Mr Steven Pryke

Consultant Orthopaedic Surgeon  
The Ipswich Hospital

### Mr Verne Johnson


Consultant Orthopaedic Surgeon  
Hull Royal Infirmary

### Mr Mark Andrews

Consultant Orthopaedic Surgeon  
Scarborough General Hospital


### Professor Alister Hart

Consultant Orthopaedic Surgeon  
Royal National Orthopaedic Hospital



## JRI Services/Education

- Research Funding
- Nurse Training
- Factory Tours
- Educational Courses
- In service support
- Consignment stock checks:
  - Implants
  - Instruments
- PACS Digital X-Ray templates



Go to [www.jri-ltd.co.uk](http://www.jri-ltd.co.uk) or  
simply scan the QR code with  
your smartphone



JRI Orthopaedics is wholly owned by Orthopaedic Research UK



**JRI Orthopaedics Ltd.**

18 Churchill Way, 35A Business Park, Chapeltown, Sheffield S35 2PY

**T:** 0114 345 0000 **F:** 0114 345 0004 **W:** [www.jri-ltd.co.uk](http://www.jri-ltd.co.uk)